

LS-78: Development of Liquid Xenon Detectors for PET

P. Amaudruz¹, F. Benard², D.A. Bryman³, C. Clements³, J. Glister¹,
L. Kurchaninov¹, J-P. Martin⁴, A. Miceli¹, F. Retiere¹, T.J. Ruth¹, V. Sossi³,
A.J. Stoessl³, H. Zhu³

¹ TRIUMF ² BC Cancer Research ³ The University of British Columbia ⁴ The University of Montreal

Detectors for PET

PET is a functional imaging technique based on detection of 511 keV annihilation photons following β^+ decay

Requirements for PET detectors:

- High sensitivity
- Sub-mm position resolution
- Good time resolution to decrease the random coincidence rate
- Good energy resolution to reject scattered photons
- High count rate capability
- Uniform response throughout the field of view
- Low cost



Existing PET systems do not satisfy all the requirements simultaneously

Applications of LXePET

Clinical applications in oncology

1. More accurately assess the effects of treatment
 - Improve patient care
 - Save costs
2. Provide better staging information
 - Increase the chance of choosing the corrected treatment

Increase sensitivity --> Reduce radioactivity dose delivered to the subject
(research and diagnosis)

Applications of LXePET

Pre-Clinical cancer studies

1. Receptor imaging: Investigation with receptor-binding tracers will be significantly enhanced

- Higher resolution --> Smaller area will be explored
- Higher sensitivity --> Impact of mass effect will be minimized

2. Heterogeneous tumors and necrosis

- Higher resolution and signal to noise ratio --> Study of heterogeneous tumors and identifying areas of central necrosis will be possible

3. Image two animals simultaneously

- Uniform resolution and high rate capability --> Imaging two animals simultaneously will have no drawbacks --> The cost of imaging studies will be reduced


LXePET

Objective: Develop a PET system that overcomes the limitations of existing PET systems and reduces detector contributions to PET to the level of intrinsic limitations

- LXePET takes advantages of the properties of liquid xenon (LXe)
- Photons interacting with the LXe produce scintillation light and ionization charge
- Simultaneous measurement of charge and light



Significant improvement to **spatial resolution, image quality, and sensitivity**



LXePET

Objective: Develop a PET system that overcomes the limitations of existing PET systems and reduces detector contributions to PET to the level of intrinsic limitations

Photon interacts with the LXe producing:


- prompt scintillation light ($\sim 2\text{ns}$)
 - measured by light detectors
- ionization
 - drifts under the applied electric field to the anode (sub-mm res. in 3D)

Anode module:

- grid of wires preceding the anode
- anode segmented into strips \perp wires
- Shielding grid




- The electron signal induced on the wires and collected by the anode strips provides the 2D position (x-y) of each interaction
- The z-coordinate is calculated from the drift time and the drift velocity



- Compton scattering reconstruction

LXePET System



12 trapezoidal sectors


Axial length = 10 cm

Drift length = 12 cm

μ-PET Sector

- TPC filled with LXe
- 32 APDs

Anode module




- Currently testing one micro-PET sector
- Designing and building new cryostat to house 2 or more sectors for coincidence PET measurements


LXePET Performance (simulated)

- MC Simulation based on Geant4, parametrization of detector response, Compton reconstruction algorithm for event reconstruction
- **Energy resolution at 511 keV = 7.5% (FWHM)**
- **Sensitivity at center of FOV = 15%**

Energy spectrum






3D sensitivity map



Small Scale Prototype Tests

- Electrons lifetime $> 100 \mu\text{s}$ (purity $< 10\text{ppb}$)
- Energy resolution 9% (FWHM)




P. Amaudruz et al., Nucl. Instr. And Meth. A (2009)

LXePET Sector Prototype


- TPC 1.1l active volume
- Anode module: 96 induction wires and 96 orthogonal anode strips
- 32 APDs
- Sector mounted inside 8.5l cryostat

2010 Tests

- Electron lifetime measurements with cosmic rays
- Coincidence measurements with Na-22 source and BrillanCe detector






Cosmic ray event – 3D visualization




Initial Measurements with Cosmic Rays

➤ Cosmic ray tracks





Electron Lifetime Measured with Cosmic Rays

- Electrons lifetime $> 70 \mu\text{s}$




Time (Z) resolution with Cosmic Ray Data

Muon track in XZplane (X = strips, Z = drift time), Strips time = Average Time




X-Position resolution from CR data





Tests with Na-22 Source

➤ 511 keV photon detected




Grid wires




Anode strip



Collimated Source



Spectrum Na²² source with Brilliance Trigger -- single strip hit events



Initial study
Very Preliminary

Non-optimum
run conditions

Sector Test: What's been accomplished so far

- Complex detector operated
- Modest purity achieved
- Position resolutions (x, z) measured:
FWHM < 1mm
- Energy resolution measured for TPC:
FWHM ~ 8% at 511 keV

Current Issues and Problems

- Purity – aim to improve by order of magnitude
 - Replace, improve liquid filters and pumps
 - Optimize liquid exchange in chamber
 - Improve cooling system for heated getter flow
- APDs
 - APD module B missing HV : fix HV connection
 - APD module A: fix 2 bad signal connections
 - APD pickup noise – investigate and fix: add coax cables in vacuum; possibly lower bandwidth
- Drift HV Connection
 - Now 10.8 KV (~ optimum for operations)
Improve -> 20 KV, desirable for tests; likely fix HV feedthrough

Schedule

Schedule

05/2010 – 12/2010

- Test prototype chamber
 - Improve purity
 - Fix APD problems
 - Study Compton reconstruction
- Complete design, manufacture, assemble cryostat, control system and two PET sectors
- Design and manufacture electronics and DAQ for multi-sector test
- Continue Geant4 simulations to optimize event reconstruction
- Develop image reconstruction algorithms

2011

- Coincidence tests with two sectors

BACKUP

LXe properties

LXe properties	
High ionization yield: 64000 e-/MeV at high E field	Large detectable ionization signal
Very small diffusion: for 1us drift 20 um	Sub-mm position resolution
High light yield: 68000 photon/MeV @ E field = 0	Bright scintillator
Fast scintillator: decay time 2.2 ns, 27 ns	Sub-ns time resolution
Wavelength of scintillation light = 178 nm	Special photo-detectors
Boiling point 165.1 K	Easy cool down
Radiation length for 511 keV photons ~ 3cm	Compact detector
Purity required < 1ppb O ₂ → e- lifetime ~ 1ms	Purification is critical

Image quality

